

REMINDERS

Specimen in Lab

When adding on lab orders to blood already drawn, please remember to change the collection priority to “Specimen in Lab” and change the collected Y/N to “yes”. If for some reason the lab cannot be added on, someone from lab will call and let you know.

These labs CANNOT be added on:

- Lactic Acid
- Ammonia Level

Merge View Feature:

Using the Merge View is a proactive duplicate check of existing orders against those that are proposed in a new PowerPlan. Turning on the Merge View allows you to view the existing orders from within your new proposed Plan. Use the Merge View before placing the PowerPlan in the Planned Status.

SURG GENERAL PREOP ADMIT (Planned Pending)

Code Status / Condition

agnosis

Use the Ctrl key to select multiple orders. Right-click and select Cancel/DC.

LIMITATIONS OF TREATMENT

Code Status / Condition (Other)

Admitting Service

Code Status

Intensity of Service (Patient Loc)

Physician Care Team

Severity of Illness (Patient Admission Category) Ordered

https://care.umh

- Renew
- Modify
- Copy
- Cancel/Reorder
- Suspend
- Activate
- Complete
- Cancel/DC
- Delete
- Add/Modify Compliance
- Print
- Disable Order Information Hyperlink

LIMITATIONS OF TREATMENT (LSE)		
Code Status / Condition (Other)		
<input type="checkbox"/>	Admitting Service	Discontinue
<input type="checkbox"/>	Code Status	Discontinue
<input type="checkbox"/>	Intensity of Service (Patient Location)	Discontinue
<input type="checkbox"/>	Physician Care Team	Discontinue
<input type="checkbox"/>	Severity of Illness (Patient Admission Category)	Discontinue

Order status is Discontinue.



These orders will be discontinued upon initiating the proposed PowerPlan.

Double Light Bulb Feature:

Using the double light bulb feature will allow you to view orders in a plan that have not yet been ordered.

Steps


1. Click on the powerplan in the plan section of the orders
 2. Click on the Double Light Bulb
- It will show all orders in the plan that have and have not been ordered.
3. You now have the option to choose orders in the plan that have not been selected (Empty box next to order). Once orders are selected, it will show as a light bulb. Finish orders, initiate and sign.



OPERATIVE REPORT AND POST OP NOTE ORDERS

ON INCORRECT ENCOUNTERS ALERT:

When patients are on SNF, LTAC, GBH, and IPR encounter locations, a new encounter called SURGOP needs to be created to enter the surgical documents and notes. This alert will fire when the notes are signed to in error the document and enter a note on the correct SURGOP encounter. If there is one created, the alert will communicate the FIN to choose.



For technical or workflow issues, please call Physician Support at 456-8142